

# Biofeedback Xpert System requirements

CF2023, June 2018

## Computer

- > PC or laptop with Pentium or compatible CPU (e.g. Celeron, Athlon), min. 2 GHz, at least 4 GB RAM
- > Graphics card (suitable for playing videos)
- > 16-bit MIDI-capable sound card with speakers or headphones
- > DVD drive
- > Hard disk (with at least 1,5 GB free space)
- > Mouse, keyboard
- > Monitor
- > A free USB port for the DONGLE
- > A free USB port for marker cable or control of an external device (if used)
- > Operating system: Windows 7, Windows 8, Windows 10

## Safety devices

If the Biofeedback is used in health care facilities the use of the following devices may be mandatory:

- > Isolating transformer for medical equipment according to EN 60601
- > Galvanic (medical) network isolation according to EN 60601 (if the computer is connected to a data network)

**Please inquire with your company's safety representative.**

Products of SCHUHFRIED GmbH are developed and in accordance with the requirements of the European Union guideline 93/42/EWG. The CE mark proves that safety-relevant regulations, EMC Standards for Medical Devices, Biocompatibility Evaluation of Medical Devices, product specific regulations and the underlying quality management system are adhered to.

**Please contact us before acquiring new equipment so that we may give you the best possible advice.**